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Standard Operating Procedure: Electrophoresis Power Supply

Device Model: PNP-1000

Manufacturer: Padideh Nozhen (Made in Iran)

1. Purpose

To provide a standardized procedure for the operation and maintenance of the **PNP-1000 Electrophoresis Power Supply** to ensure accurate separation of nucleic acids.

2. Scope of Operation

Electrophoresis is primarily used for the separation of **DNA and RNA** fragments. The principle of separation relies on the migration of charged molecules through a matrix (agarose gel) under the influence of an electric field. The molecules move from the **negative electrode (cathode)** toward the **positive electrode (anode)** within a specific running buffer contained in the electrophoresis tank.

3. Responsibility

The responsibility for executing this protocol lies with the relevant laboratory experts and technicians.

4. Operating Procedure

A. Gel Preparation

- Agarose Casting:** Dissolve agarose powder in the required buffer by boiling until completely clear.
- Sealing:** Seal the gel tray/casting tray using waterproof tape to prevent leakage.
- Comb Placement:** Position the comb so that it sits approximately **1 mm** above the bottom of the tray, then pour the molten agarose into the tray.
- Comb Removal:** Once the gel has solidified, carefully remove the comb and the sealing tape from around the gel chamber.

B. Sample Loading and Setup

- Buffer Addition:** Fill the electrophoresis tank with **TBE buffer (0.5X - 1X)** until the buffer level is approximately **0.5 cm** above the surface of the gel.
- Loading Samples:** Load the samples into the gel wells along with an appropriate **molecular weight marker (DNA Ladder)**.
- Note on Loading Dye:** If the PCR product does not contain a pre-mixed dye, use a **6X Loading Dye** (which contains DNA density agents). Mix the sample and dye at a ratio of **5 parts PCR product to 1 part Loading Dye**.

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- Power On & Pre-setting:** Turn on the power switch located at the rear of the device. **Before connecting the cables to the tank**, adjust the voltage settings to between **90V and 110V**.
- Tank Connection:** Place the gel inside the tank and add the running buffer (ensuring it covers the gel by 0.5 cm). Close the tank lid securely.
- Electrode Connection:** Connect the electrodes to the tank.
 - Black Cable** → \rightarrow → **Negative Pole (Cathode)**
 - Red Cable** → \rightarrow → **Positive Pole (Anode)**
- Crucial Directional Note:** Because DNA is negatively charged, it migrates from the **Cathode (-)** toward the **Anode (+)**. Therefore, the wells must be positioned at the **Black (negative)** end of the tank to ensure the DNA runs toward the Red (positive) end.

D. Post-Run Procedures

- Power Down:** Once the electrophoresis is complete, turn off the power supply first. Disconnect the cables by **grasping the plug heads** rather than pulling the wires.
- Gel Retrieval:** Remove the tank lid and carefully transfer the gel to an imaging system (Gel Doc) for visualization.
- Cleaning:** Empty the buffer from the tank. Thoroughly rinse the tray, comb, and tank with **distilled water**.

Warning: Failure to promptly empty and clean the electrophoresis tank, or the continued use of old, precipitate-heavy buffers, will cause the internal connection wires to lose their conductivity over time.

5. Safety Instructions

- PPE:** Always wear appropriate laboratory gloves during operation.
- Stability:** Never attempt to move or reposition the power supply unit while it is switched ON.
- Voltage Control:** Avoid connecting the tank to voltages exceeding the recommended limits. Excessive voltage causes buffer heating, which degrades both the gel integrity and the sample quality.
- Electrical Hazard:** To prevent electric shock, **never** attempt to add or remove a gel from the tank while the power supply is connected and turned on.
- Lid Safety:** Always keep the tank lid closed during operation to prevent accidental contact with the electric field.
- Liquid Spills:** Prevent any liquid from spilling onto the power supply unit, as this can damage the electronic circuit boards and significantly increase the risk of electric shock.
- Contamination Control:** To minimize cross-contamination, do not move electrophoresis equipment and accessories to other departments within the laboratory.

6. Maintenance

- Conductivity Preservation:** As noted previously, timely drainage and cleaning of the tank are mandatory. Using aged or sediment-heavy buffers leads to the degradation of the tank's internal wiring and reduces conductivity.

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2. **Cable Handling:** To prevent wire breakage or electrical interruptions, always grasp the **connector head** when plugging or unplugging cables from the tank.
3. **Drying:** At the end of every work session, rinse all utilized components (tank, tray, comb) with distilled water and ensure they are completely dried before storage.

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